

For Research Use Only

EGR3 Polyclonal antibody, PBS Only

Catalog Number: 26495-1-PBS



Basic Information

Catalog Number:

26495-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24152

GenBank Accession Number:

BC104765

GeneID (NCBI):

1960

UNIPROT ID:

Q06889

Full Name:

early growth response 3

Calculated MW:

387 aa, 43 kDa

Observed MW:

43 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

The early growth response (Egr) family of zinc finger transcription factors regulates critical genetic programs in cellular growth, differentiation, and function. EGR3 (Early Growth Response 3; also known as PILOT) is a member of the early growth response (Egr) gene family of transcription factors. EGR3 regulates the expression of genes involved in inflammation, including genes encoding various secreted cytokines and protease inhibitors. EGR3 functions as a master regulator of genes differentially expressed in the brains of patients with neuropsychiatric illnesses ranging from schizophrenia and bipolar disorder to Alzheimer's disease.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26495-1-AP (EGR3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26495-1-PBS in a different storage buffer formulation.