

For Research Use Only

ELOVL7 Polyclonal antibody

Catalog Number: 26615-1-AP



Basic Information

Catalog Number: 26615-1-AP	GenBank Accession Number: BC130310	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 800 µg/ml by Nanodrop;	GeneID (NCBI): 79993	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: ELOVL family member 7, elongation of long chain fatty acids (yeast)	
Isotype: IgG	Calculated MW: 281 aa, 33 kDa	
Immunogen Catalog Number: AG25008	Observed MW: 33 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells, HepG2 cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26615-1-AP (ELOVL7 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.