

For Research Use Only

WSTF Polyclonal antibody

Catalog Number: 26641-1-AP



Basic Information

Catalog Number:

26641-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24403

GenBank Accession Number:

BC136520

GeneID (NCBI):

9031

UNIPROT ID:

Q9UIGO

Full Name:

bromodomain adjacent to zinc finger domain, 1B

Calculated MW:

1483 aa, 171 kDa

Observed MW:

170 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : SH-SY5Y cells, HEK-293 cells, HeLa cells, HepG2 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26641-1-AP (WSTF antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.