

For Research Use Only

CACNA2D4 Polyclonal antibody, PBS Only

Catalog Number: 26674-1-PBS



Basic Information

Catalog Number:

26674-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24464

GenBank Accession Number:

BC150186

GeneID (NCBI):

93589

UNIPROT ID:

Q7Z3S7

Full Name:

calcium channel, voltage-dependent, alpha 2/delta subunit 4

Calculated MW:

1137 aa, 128 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

Calcium voltage-gated channel auxiliary subunit alpha2delta 4 (CACNA2D4, also known as RCD4) is a member of the alpha-2/delta subunit family within the voltage-dependent calcium channel complex. The $\alpha 2\delta 4$ subunit of retinal voltage-gated calcium channels (Cav), is crucial for developing and mature exocytotic function of the photoreceptor ribbon synapse (PMID: 29875267). CACNA2D4 is localized to the retina, particularly in photoreceptors, where it is essential for organizing synaptic ribbons and setting Cav1.4 voltage sensitivity (PMID: 28262416). In humans, mutations in CACNA2D4 are associated with autosomal recessive cone dystrophy, a condition that affects the photoreceptor cells in the retina (PMID: 16877424).

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26674-1-AP (CACNA2D4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26674-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 26674-1-AP (CACNA2D4 antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26674-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 26674-1-AP (CACNA2D4 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 26674-1-PBS in a different storage buffer formulation.