

For Research Use Only

# ABCD3 Polyclonal antibody

Catalog Number: 26720-1-AP



## Basic Information

<b>Catalog Number:</b> 26720-1-AP	<b>GenBank Accession Number:</b> BC068509	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 150ul , Concentration: 700 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 5825	<b>Recommended Dilutions:</b> WB 1:1000-1:4000
<b>Source:</b> Rabbit	<b>Full Name:</b> ATP-binding cassette, sub-family D (ALD), member 3	
<b>Isotype:</b> IgG	<b>Observed MW:</b> 70 kDa	
<b>Immunogen Catalog Number:</b> AG24760		

## Applications

**Tested Applications:**  
WB, ELISA

**Species Specificity:**  
human, mouse, rat

### Positive Controls:

WB : mouse colon tissue, mouse kidney tissue, mouse liver tissue, rat liver tissue

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com      W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26720-1-AP (ABCD3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.