

For Research Use Only

ABCD3 Polyclonal antibody

Catalog Number: 26720-1-AP **1 Publications**



Basic Information

Catalog Number: 26720-1-AP	GenBank Accession Number: BC068509	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 700 ug/ml by Nanodrop;	GeneID (NCBI): 5825	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: P28288	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family D (ALD), member 3	
Immunogen Catalog Number: AG24760	Observed MW: 70 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls:
Cited Applications: WB	WB : mouse colon tissue, mouse kidney tissue, mouse liver tissue, rat liver tissue
Species Specificity: human, mouse, rat	
Cited Species: human	

Notable Publications

Author	Pubmed ID	Journal	Application
Meng Liu	38649521	Cell Mol Life Sci	WB

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26720-1-AP (ABCD3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.