

For Research Use Only

PLB1 Polyclonal antibody, PBS Only

Catalog Number: 26835-1-PBS



Basic Information

Catalog Number:

26835-1-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25240

GenBank Accession Number:

BC150633

GeneID (NCBI):

151056

UNIPROT ID:

Q6P1J6

Full Name:

phospholipase B1

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, IHC, ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



DU 145 cells were subjected to SDS PAGE followed by western blot with 26835-1-AP (PLB1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26835-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate hyperplasia tissue slide using 26835-1-AP (PLB1 Antibody) at dilution of 1:200 (under 10x lens). This data was developed using the same antibody clone with 26835-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human prostate hyperplasia tissue slide using 26835-1-AP (PLB1 Antibody) at dilution of 1:200 (under 40x lens). This data was developed using the same antibody clone with 26835-1-PBS in a different storage buffer formulation.