

For Research Use Only

SLC4A3 Polyclonal antibody

Catalog Number: 26883-1-AP

1 Publications



Basic Information

Catalog Number:

26883-1-AP

Size:

150ul, Concentration: 1000 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25457

GenBank Accession Number:

BC146656

GeneID (NCBI):

6508

UNIPROT ID:

P48751

Full Name:

solute carrier family 4, anion exchanger, member 3

Calculated MW:

1232 aa, 136 kDa

Observed MW:

135-140 kDa, 104 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse brain tissue, mouse heart tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Eunus S Ali	35772404	Mol Cell	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26883-1-AP (SLC4A3 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.