

For Research Use Only

# HNRNPM Polyclonal antibody

Catalog Number: 26897-1-AP

Featured Product

8 Publications



## Basic Information

Catalog Number:

26897-1-AP

Size:

150ul, Concentration: 400 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25474

GenBank Accession Number:

BC000138

GeneID (NCBI):

4670

UNIPROT ID:

P52272

Full Name:

heterogeneous nuclear ribonucleoprotein M

Calculated MW:

78 kDa

Observed MW:

73-77 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

IHC 1:200-1:800

IF/ICC 1:50-1:500

## Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB, IF, IP

Species Specificity:

human

Cited Species:

human, mouse

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

Positive Controls:

WB : COLO 320 cells, MCF-7 cells, HeLa cells, NIH/3T3 cells, Neuro-2a cells

IHC : human tonsillitis tissue, human colon

IF/ICC : HepG2 cells,

## Background Information

HNRNPM, also named as NAGR1, is a 730 amino acid protein, which localizes in nucleus. HNRNPM as a pre-mRNA binding protein in vivo, binds avidly to poly(G) and poly(U) RNA homopolymers in vitro. HNRNPM is involved in splicing. HNRNPM acts as a receptor for carcinoembryonic antigen in Kupffer cells, may initiate a series of signaling events leading to tyrosine phosphorylation of proteins and induction of IL-1 alpha, IL-6, IL-10 and tumor necrosis factor alpha cytokines. HNRNPM exists two isoforms with MV 77 kDa and 73 kDa.

## Notable Publications

Author	Pubmed ID	Journal	Application
Tao Wang	36151083	Nat Commun	IF
Wen-Juan Li	33782401	Nat Commun	WB
Pingfu Hou	31291975	J Exp Clin Cancer Res	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Selected Validation Data




COLO 320 cells were subjected to SDS PAGE followed by western blot with 26897-1-AP (HNRNPM antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.




WB result of HNRNPM antibody (26897-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-HNRNPM transfected MCF-7 cells.




COLO 320 cells were subjected to SDS PAGE followed by western blot with 26897-1-AP (HNRNPM antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human tonsillitis tissue slide using 26897-1-AP (HNRNPM antibody) at dilution of 1:400 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using HNRNPM antibody (26897-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L), CL594-Phalloidin (red).



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using HNRNPM antibody (26897-1-AP) at dilution of 1:400 and CoraLite®594-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-4), CL488-Phalloidin (green).