

For Research Use Only

RIT1 Polyclonal antibody, PBS Only

Catalog Number: 26977-1-PBS



Basic Information

Catalog Number: 26977-1-PBS	GenBank Accession Number: BC104186	Purification Method: Antigen affinity purification
Size: 100ug, Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6016	
Source: Rabbit	UNIPROT ID: Q92963	
Isotype: IgG	Full Name: Ras-like without CAAX 1	
Immunogen Catalog Number: AG25636	Calculated MW: 25 kDa	
	Observed MW: 25-28 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

RIT1 (Ras-like without CAAX protein 1) is also named as RIBB, ROC1 and GTP-binding protein Rit1. RIT1 serves as a regulatory factor for neuronal cell proliferation, survival, and differentiation (PMID: 28007959). Rit1 mediates oxidative stress resistance, contributing to cell survival via the p38 MAPK signaling pathway (PMID: 21737674). As a small G protein of Ras family, RIT1 plays a critical role in various tumors (such as hepatocellular carcinoma, HCC) (PMID: 38017479). It is worth noting that RIT1 is frequently amplified in various human cancers including HCC, lung adenocarcinoma, cholangiocarcinoma, uterine carcinosarcoma, breast cancer, and ovarian cancer (PMID: 38017479).

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 26977-1-AP (RIT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 26977-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Methanol) fixed MCF-7 cells using RIT1 antibody (26977-1-AP) at dilution of 1:200 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 26977-1-PBS in a different storage buffer formulation.