

For Research Use Only

TAF5 Polyclonal antibody

Catalog Number: 27041-1-AP



Basic Information

Catalog Number:

27041-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25742

GenBank Accession Number:

BC136348

GeneID (NCBI):

6877

UNIPROT ID:

Q15542

Full Name:

TAF5 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 100kDa

Calculated MW:

87 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HCT 116 cells, HEK-293 cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27041-1-AP (TAF5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.