

For Research Use Only

# TAF5 Polyclonal antibody, PBS Only

Catalog Number: 27041-1-PBS



## Basic Information

Catalog Number:

27041-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25742

GenBank Accession Number:

BC136348

GeneID (NCBI):

6877

UNIPROT ID:

Q15542

Full Name:

TAF5 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 100kDa

Calculated MW:

87 kDa

Observed MW:

100 kDa

Purification Method:

Antigen affinity purification

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27041-1-AP (TAF5 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27041-1-PBS in a different storage buffer formulation.



WB result of TAF5 antibody (27041-1-AP; 1:600; incubated at room temperature for 1.5 hours) with sh-Control and sh-TAF5 transfected HEK-293 cells. This data was developed using the same antibody clone with 27041-1-PBS in a different storage buffer formulation.