

For Research Use Only

APOC2 Polyclonal antibody

Catalog Number: 27045-1-AP



Basic Information

Catalog Number:

27045-1-AP

Size:

150ul , Concentration: 1000 ug/ml by Nanodrop and 667 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24650

GenBank Accession Number:

BC005348

GeneID (NCBI):

344

UNIPROT ID:

P02655

Full Name:

apolipoprotein C-II

Observed MW:

9 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : human plasma tissue,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human plasma were subjected to SDS PAGE followed by western blot with 27045-1-AP (APOC2 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.