

For Research Use Only

ZRANB1 Polyclonal antibody, PBS Only

Catalog Number: 27050-1-PBS



Basic Information

Catalog Number:

27050-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25789

GenBank Accession Number:

AL832925

GeneID (NCBI):

54764

UNIPROT ID:

Q9UGI0

Full Name:

zinc finger, RAN-binding domain containing 1

Calculated MW:

81 kDa

Observed MW:

81 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

IP, Indirect ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



IP result of anti-ZRANB1 (IP:27050-1-AP, 4ug; Detection:27050-1-AP 1:500) with HeLa cells lysate 1600 ug. This data was developed using the same antibody clone with 27050-1-PBS in a different storage buffer formulation.



IP result of anti-ZRANB1 (IP:27050-1-AP, 4ug; Detection:27050-1-AP 1:500) with mouse cerebellum tissue lysate 2000 ug. This data was developed using the same antibody clone with 27050-1-PBS in a different storage buffer formulation.