

For Research Use Only

Cytokeratin 73 Polyclonal antibody, PBS Only

Catalog Number: 27055-1-PBS



Basic Information

Catalog Number: 27055-1-PBS	GenBank Accession Number: BC109212	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 319101	
Source: Rabbit	UNIPROT ID: Q86Y46	
Isotype: IgG	Full Name: keratin 73	
Immunogen Catalog Number: AG25716	Observed MW: 59 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, rat

Background Information

Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. This gene encodes a protein that is expressed in the inner root sheath of hair follicles. The type II keratins are clustered in a region of chromosome 12q13.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



rat skin tissue were subjected to SDS PAGE followed by western blot with 27055-1-AP (KRT73 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27055-1-PBS in a different storage buffer formulation.