

For Research Use Only

PSRC1 Polyclonal antibody

Catalog Number: 27130-1-AP



Basic Information

Catalog Number: 27130-1-AP	GenBank Accession Number: BC056909	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 900 µg/ml by Nanodrop and 467 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 84722	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: proline/serine-rich coiled-coil 1	
Isotype: IgG	Calculated MW: 39 kDa	
Immunogen Catalog Number: AG26039	Observed MW: 39 kDa	

Applications

Tested Applications: WB,ELISA	Positive Controls: WB : A431 cells, HEK-293T cells
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A431 cells were subjected to SDS PAGE followed by western blot with 27130-1-AP (PSRC1 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.