

For Research Use Only

# ABCB8 Polyclonal antibody

Catalog Number: 27158-1-AP



## Basic Information

|                                      |  |   |
|--------------------------------------|--|---|
| Catalog Number:<br>27158-1-AP        | GenBank Accession Number:<br>BC141814                                | Purification Method:<br>Antigen affinity purification |
| Size:<br>150ul                       | GeneID (NCBI):<br>11194  |   |
| Source:<br>Rabbit                    | Full Name:<br>ATP-binding cassette, sub-family B (MDR/TAP), member 8 |   |
| Isotype:<br>IgG                      | Calculated MW:<br>735 aa, 80 kDa                                     |   |
| Immunogen Catalog Number:<br>AG25935 |  |   |

## Applications

Tested Applications:  
ELISA

Species Specificity:  
human, mouse

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data