

For Research Use Only

LMBRD1 Polyclonal antibody, PBS Only

Catalog Number: 27203-1-PBS



Basic Information

Catalog Number:

27203-1-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24478

GenBank Accession Number:

BC010360

GeneID (NCBI):

55788

UNIPROT ID:

Q9NUN5

Full Name:

LMBR1 domain containing 1

Calculated MW:

540 aa, 61 kDa

Observed MW:

47 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with 27203-1-AP (LMBRD1 antibody) at dilution of 1:400 incubated at room temperature for 1.5 hours.



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with 27203-1-AP (LMBRD1 Antibody) at dilution of 1:100 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27203-1-PBS in a different storage buffer formulation.