

For Research Use Only

RAB5C Polyclonal antibody

Catalog Number: 27219-1-AP

1 Publications



Basic Information

Catalog Number:

27219-1-AP

Size:

150ul, Concentration: 900 ug/ml by Nanodrop and 467 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25849

GenBank Accession Number:

BC106039

GeneID (NCBI):

5878

UNIPROT ID:

P51148

Full Name:

RAB5C, member RAS oncogene family

Calculated MW:

216 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

human

Positive Controls:

WB : A549 cells, U-937 cells

Notable Publications

Author	Pubmed ID	Journal	Application
Jun Yan	34367131	Front Immunol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



A549 cells were subjected to SDS PAGE followed by western blot with 27219-1-AP (RAB5C antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.