

For Research Use Only

PAX2 Polyclonal antibody

Catalog Number: 27222-1-AP



Basic Information

| | | |
|---------------------------------------------|----------------------------------------------|--------------------------------------------------------------|
| Catalog Number: 27222-1-AP | GenBank Accession Number: BC148710 | Purification Method: Antigen affinity purification |
| Size: 150ul | GeneID (NCBI): 5076 | Recommended Dilutions: IHC 1:50-1:500 |
| Source: Rabbit | Full Name: paired box 2 | |
| Isotype: IgG | Calculated MW: 431 aa, 46 kDa | |
| Immunogen Catalog Number: AG25929 | | |

Applications

| | |
|------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
| Tested Applications: IHC, ELISA | Positive Controls: IHC : human renal cell carcinoma tissue, |
| Species Specificity: human | |
| Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0 | |

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data