

For Research Use Only

KANSL2 Polyclonal antibody

Catalog Number: 27261-1-AP



Basic Information

Catalog Number: 27261-1-AP	GenBank Accession Number: BC009746	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 850 µg/ml by Nanodrop and 400 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 54934	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: chromosome 12 open reading frame 41	
Isotype: IgG	Observed MW: 50 kDa~55 kDa	
Immunogen Catalog Number: AG25983		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse heart tissue,
Species Specificity: human, mouse	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 27261-1-AP (KANSL2) antibody at dilution of 1:2000 incubated at room temperature for 1.5 hours.