

For Research Use Only

BARX1 Polyclonal antibody

Catalog Number: 27281-1-AP **1 Publications**



Basic Information

Catalog Number: 27281-1-AP	GenBank Accession Number: BC064363	Purification Method: Antigen affinity purification
Size: 150ul	GeneID (NCBI): 56033	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	Full Name: BARX homeobox 1	
Isotype: IgG		
Immunogen Catalog Number: AG25833		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse heart tissue,
Cited Applications: WB	
Species Specificity: human	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Danjie Liu	35859853	Dose Response	WB

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 27281-1-AP (BARX1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.