

For Research Use Only

BARX1 Polyclonal antibody

Catalog Number: 27281-1-AP

1 Publications



Basic Information

Catalog Number:

27281-1-AP

Size:

150ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25833

GenBank Accession Number:

BC064363

GeneID (NCBI):

56033

UNIPROT ID:

Q9HBU1

Full Name:

BARX homeobox 1

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human

Cited Species:

mouse

Positive Controls:

WB : mouse heart tissue,

Notable Publications

Author	Pubmed ID	Journal	Application
Danjie Liu	35859853	Dose Response	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 27281-1-AP (BARX1 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.