

For Research Use Only

C14orf153 Polyclonal antibody

Catalog Number: 27300-1-AP **1 Publications**



Basic Information

| | | |
|--|---|--|
| Catalog Number: 27300-1-AP | GenBank Accession Number: BC007412 | Purification Method: Antigen affinity purification |
| Size: 150ul , Concentration: 267 µg/ml by Bradford method using BSA as the standard; | GeneID (NCBI): 84334 | |
| Source: Rabbit | Full Name: chromosome 14 open reading frame 153 | |
| Isotype: IgG | | |
| Immunogen Catalog Number: AG23870 | | |

Applications

Tested Applications:
IF, ELISA

Cited Applications:
IF, WB

Species Specificity:
human

Cited Species:
human, mouse

Background Information

Notable Publications

| Author | Pubmed ID | Journal | Application |
|-------------|-----------|--------------|-------------|
| Alba Signes | 30552096 | EMBO Mol Med | IF, WB |

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data