

For Research Use Only

DNAJC5G Polyclonal antibody

Catalog Number: 27373-1-AP



Basic Information

Catalog Number:

27373-1-AP

Size:

150ul, Concentration: 570 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11748

GenBank Accession Number:

BC035194

GeneID (NCBI):

285126

UNIPROT ID:

Q8N752

Full Name:

DnaJ (Hsp40) homolog, subfamily C, member 5 gamma

Calculated MW:

189 aa, 21 kDa

Observed MW:

21 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PC-3 cells, SH-SY5Y cells, SK-N-SH cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27373-1-AP (DNAJC5G antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.