

For Research Use Only

CNTNAP3 Polyclonal antibody

Catalog Number: 27410-1-AP



Basic Information

Catalog Number:
27410-1-AP

Size:
150ul , Concentration: 500 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;

Source:
Rabbit

Isotype:
IgG

Immunogen Catalog Number:
AG26624

GenBank Accession Number:
BC132737

GeneID (NCBI):
79937

Full Name:
contactin associated protein-like 3

Observed MW:
124 kDa

Purification Method:
Antigen affinity purification

Recommended Dilutions:
WB 1:500-1:2000

Applications

Tested Applications:
WB, ELISA

Species Specificity:
human

Positive Controls:

WB : HEK-293 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

***** 20ul sizes contain 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 27410-1-AP (CNTNAP3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.