

For Research Use Only

WDR91 Polyclonal antibody

Catalog Number: 27424-1-AP



Basic Information

Catalog Number:

27424-1-AP

Size:

150ul , Concentration: 750 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26697

GenBank Accession Number:

BC017246

GeneID (NCBI):

29062

UNIPROT ID:

A4D1P6

Full Name:

WD repeat domain 91

Observed MW:

42 kDa, 83 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, Jurkat cells, SH-SY5Y cells, U2OS cells

Background Information

WDR91 (1-405 aa) can be detected 42 kDa (PMID:28860274).

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27424-1-AP (WDR91 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.