

For Research Use Only

WDR91 Polyclonal antibody

Catalog Number: 27424-1-AP



Basic Information

Catalog Number: 27424-1-AP	GenBank Accession Number: BC017246	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 750 µg/ml by Nanodrop;	GeneID (NCBI): 29062	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: WD repeat domain 91	
Isotype: IgG	Observed MW: 42 kDa, 83 kDa	
Immunogen Catalog Number: AG26697		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, Jurkat cells, SH-SY5Y cells, U2OS cells
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27424-1-AP (WDR91 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.