

For Research Use Only

SH3BP5L Polyclonal antibody, PBS Only

Catalog Number: 27446-1-PBS



Basic Information

Catalog Number:

27446-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23990

GenBank Accession Number:

BC017254

GeneID (NCBI):

80851

UNIPROT ID:

Q7L8J4

Full Name:

SH3-binding domain protein 5-like

Observed MW:

65-70 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

SH3BP5L (SH3 Binding Domain Protein 5 Like) is a protein-coding gene that encodes a protein with a conserved N-terminus and functions as a guanine nucleotide exchange factor (GEF) for Rab11 GTPases.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27446-1-AP (SH3BP5L antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27446-1-PBS in a different storage buffer formulation.