

For Research Use Only

C21orf70 Polyclonal antibody

Catalog Number: 27447-1-AP



Basic Information

Catalog Number: 27447-1-AP	GenBank Accession Number: BC009341	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop and 267 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 85395	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: chromosome 21 open reading frame 70	
Isotype: IgG	Observed MW: 24 kDa	
Immunogen Catalog Number: AG24033		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HeLa cells,
Species Specificity: human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 27447-1-AP (C21orf70 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.