

For Research Use Only

PEX26 Polyclonal antibody

Catalog Number: 27472-1-AP **2 Publications**



Basic Information

Catalog Number: 27472-1-AP	GenBank Accession Number: BC016280	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 55670	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	Full Name: peroxisomal biogenesis factor 26	
Isotype: IgG	Calculated MW: 34 kDa	
Immunogen Catalog Number: AG26529	Observed MW: 30-34 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A431 cells, HEK-293 cells
Cited Applications: IHC, WB	
Species Specificity: human	
Cited Species: human	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Richard G Lee	38129691	Nat Cell Biol	WB
Bianbian Yan	37957408	Dig Dis Sci	IHC, WB

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27472-1-AP (PEX26 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.