

For Research Use Only

ADAM22 Polyclonal antibody, PBS Only

Catalog Number: 27496-1-PBS



Basic Information

Catalog Number:

27496-1-PBS

Size:

100ug, Concentration: 1 mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26783

GenBank Accession Number:

BC125118

GeneID (NCBI):

53616

UNIPROT ID:

Q9P0K1

Full Name:

ADAM metalloproteinase domain 22

Observed MW:

75-90 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27496-1-AP (ADAM22 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27496-1-PBS in a different storage buffer formulation.