

For Research Use Only

HES1 Polyclonal antibody

Catalog Number: 27542-1-AP

4 Publications



Basic Information

Catalog Number:

27542-1-AP

Size:

150ul, Concentration: 367 µg/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26393

GenBank Accession Number:

BC039152

GeneID (NCBI):

3280

UNIPROT ID:

Q14469

Full Name:

hairy and enhancer of split 1, (Drosophila)

Calculated MW:

30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human, mouse, rat

Notable Publications

Author	Pubmed ID	Journal	Application
Tianrong Zhang	31391046	Lipids Health Dis	WB
Xueying Jin	38703204	Cell Mol Life Sci	WB
Menghuan Wang	38134973	J Nutr Biochem	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

