

For Research Use Only

TIAM1 Polyclonal antibody

Catalog Number: 27694-1-AP

1 Publications



Basic Information

Catalog Number:

27694-1-AP

Size:

150ul, Concentration: 1200 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26777

GenBank Accession Number:

NM_003253

GeneID (NCBI):

7074

UNIPROT ID:

Q13009

Full Name:

T-cell lymphoma invasion and metastasis 1

Calculated MW:

178 kDa

Observed MW:

178 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:16000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, mouse

Cited Species:

human

Positive Controls:

WB : HEK-293T cells, HCT 116 cells, Neuro-2a cells

Notable Publications

Author	Pubmed ID	Journal	Application
Tianyu Han	38460155	Adv Sci (Weinh)	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27694-1-AP (TIAM1 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.