

For Research Use Only

# FAM91A1 Polyclonal antibody

Catalog Number: 27738-1-AP

1 Publications



## Basic Information

Catalog Number:

27738-1-AP

Size:

150ul, Concentration: 350 ug/ml by Nanodrop and 267 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26902

GenBank Accession Number:

BC030520

GeneID (NCBI):

157769

UNIPROT ID:

Q658Y4

Full Name:

family with sequence similarity 91, member A1

Observed MW:

94 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

IP

Species Specificity:

Human

Cited Species:

human

Positive Controls:

WB : A549 cells,

## Notable Publications

| Author       | Pubmed ID | Journal | Application |
|--------------|-----------|---------|-------------|
| Huaqing Deng | 39013469  | Cell    | IP          |

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



A549 and HeLa cells were subjected to SDS PAGE followed by western blot with 27738-1-AP (FAM91A1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.