

For Research Use Only

HAPLN2 Polyclonal antibody

Catalog Number: 27818-1-AP



Basic Information

Catalog Number: 27818-1-AP	GenBank Accession Number: BC029864	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 350 µg/ml by Nanodrop and 267 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 60484	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	Full Name: hyaluronan and proteoglycan link protein 2	
Isotype: IgG	Calculated MW: 38 kDa	
Immunogen Catalog Number: AG27169	Observed MW: 35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : pig brain tissue, rat brain tissue
Species Specificity: Human, Rat, Pig	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig brain tissue were subjected to SDS PAGE followed by western blot with 27818-1-AP (HAPLN2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.