

For Research Use Only

HAPLN2 Polyclonal antibody

Catalog Number: 27818-1-AP



Basic Information

Catalog Number:

27818-1-AP

Size:

150ul , Concentration: 350 ug/ml by Nanodrop and 267 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27169

GenBank Accession Number:

BC029864

GeneID (NCBI):

60484

UNIPROT ID:

Q9GZV7

Full Name:

hyaluronan and proteoglycan link protein 2

Calculated MW:

38 kDa

Observed MW:

35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Rat, Pig

Positive Controls:

WB : pig brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



pig brain tissue were subjected to SDS PAGE followed by western blot with 27818-1-AP (HAPLN2 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.