

For Research Use Only

RPUSD4 Polyclonal antibody

Catalog Number: 27843-1-AP



Basic Information

Catalog Number:

27843-1-AP

Size:

150ul , Concentration: 400 ug/ml by Nanodrop and 333 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27225

GenBank Accession Number:

BC014131

GeneID (NCBI):

84881

UNIPROT ID:

Q96CM3

Full Name:

RNA pseudouridylate synthase domain containing 4

Observed MW:

39-42 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse, Rat

Positive Controls:

WB : mouse liver tissue, rat liver tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27843-1-AP (RPUSD4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.