

For Research Use Only

RyDEN Polyclonal antibody

Catalog Number: 27865-1-AP

Featured Product

2 Publications



Basic Information

Catalog Number:

27865-1-AP

Size:

150ul, Concentration: 500 ug/ml by Nanodrop and 267 ug/ml by Bradford method using BSA as the standard;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26838

GenBank Accession Number:

BC035817

GeneID (NCBI):

55337

UNIPROT ID:

Q9NUL5

Full Name:

chromosome 19 open reading frame

66

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, Mouse, Rat

Cited Species:

human, mouse

Positive Controls:

WB: THP-1 cells, mouse liver tissue, rat liver tissue

Notable Publications

Author	Pubmed ID	Journal	Application
Matthias M Zimmer	34893599	Nat Commun	WB
Natasha W Hanners	34873063	Proc Natl Acad Sci U S A	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



THP-1 cells were subjected to SDS PAGE followed by western blot with 27865-1-AP (RyDEN antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 27865-1-AP (RyDEN antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.