

For Research Use Only

UBE2J1 Polyclonal antibody, PBS Only

Catalog Number:27890-1-PBS



Basic Information

Catalog Number:

27890-1-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27310

GenBank Accession Number:

BC013973

GeneID (NCBI):

51465

UNIPROT ID:

Q9Y385

Full Name:

ubiquitin-conjugating enzyme E2, J1 (UBC6 homolog, yeast)

Calculated MW:

318 aa, 35 kDa

Observed MW:

35-40 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 27890-1-AP (UBE2J1 antibody) at dilution of 1:30000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27890-1-PBS in a different storage buffer formulation.