

For Research Use Only

C3orf15 Polyclonal antibody

Catalog Number: 27967-1-AP



Basic Information

Catalog Number:

27967-1-AP

Size:

150ul, Concentration: 550 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27420

GenBank Accession Number:

BC126394

GeneID (NCBI):

89876

UNIPROT ID:

Q7Z4T9

Full Name:

chromosome 3 open reading frame 15

Observed MW:

11 kDa, 71 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB: mouse liver tissue, mouse testis tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27967-1-AP (C3orf15 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.