

For Research Use Only

C3orf15 Polyclonal antibody

Catalog Number: 27967-1-AP



Basic Information

Catalog Number: 27967-1-AP	GenBank Accession Number: BC126394	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 550 µg/ml by Nanodrop;	GeneID (NCBI): 89876	Recommended Dilutions: WB 1:500-1:3000
Source: Rabbit	Full Name: chromosome 3 open reading frame 15	
Isotype: IgG	Observed MW: 11 kDa, 71 kDa	
Immunogen Catalog Number: AG27420		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse liver tissue, mouse testis tissue
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27967-1-AP (C3orf15 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.