

For Research Use Only

RBMXL2 Polyclonal antibody, PBS Only

Catalog Number: 27987-1-PBS



Basic Information

Catalog Number: 27987-1-PBS	GenBank Accession Number: BC057796	Purification Method: Antigen affinity purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 27288	
Source: Rabbit	UNIPROT ID: O75526	
Isotype: IgG	Full Name: RNA binding motif protein, X-linked-like 2	
Immunogen Catalog Number: AG27624	Calculated MW: 43 kDa	
	Observed MW: 43 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

RBMXL2 (RNA-binding motif protein, X-linked-like-2), also known as heterogeneous nuclear ribonucleoproteins G-testis or hnRNP-GT, is a testis-specific nuclear RNA binding protein that plays a crucial role in male meiosis (PMID: 34927545). The protein is part of an anciently diverged family of RNA-binding proteins and shares significant homology with RBMX, another RNA-binding protein that is important for controlling splicing, transcription, and genome stability (PMID: 39356106). The expression of RBMXL2 is restricted to the testis, and it is highly expressed during the pachytene and diplotene stages of meiosis, which are characterized by high levels of transcription and alternative splicing (PMID: 34927545).

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27987-1-AP (RBMXL2 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 27987-1-PBS in a different storage buffer formulation.