

For Research Use Only

SCN1A Polyclonal antibody, PBS Only

Catalog Number: 28079-1-PBS



Basic Information

Catalog Number:

28079-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25942

GenBank Accession Number:

NM_001165963

GeneID (NCBI):

6323

UNIPROT ID:

P35498

Full Name:

sodium channel, voltage-gated, type I, alpha subunit

Calculated MW:

229 kDa

Observed MW:

229-240 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PC-12 cells were subjected to SDS PAGE followed by western blot with 28079-1-AP (SCN1A antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 28079-1-PBS in a different storage buffer formulation.