

For Research Use Only

DLGAP1 Polyclonal antibody

Catalog Number: 28124-1-AP



Basic Information

Catalog Number:

28124-1-AP

Size:

150ul , Concentration: 450 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27182

GenBank Accession Number:

BC070175

GeneID (NCBI):

9229

UNIPROT ID:

O14490

Full Name:

discs, large (Drosophila) homolog-associated protein 1

Calculated MW:

543 aa, 60 kDa

Observed MW:

110-130 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : mouse brain tissue, rat brain tissue, rat cerebellum tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28124-1-AP (DLGAP1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.