

For Research Use Only

APH1B Polyclonal antibody

Catalog Number: 28175-1-AP



Basic Information

Catalog Number:

28175-1-AP

Size:

150ul , Concentration: 300 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14066

GenBank Accession Number:

BC020905

GeneID (NCBI):

83464

UNIPROT ID:

Q8WW43

Full Name:

anterior pharynx defective 1 homolog B (C. elegans)

Calculated MW:

257 aa, 28 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : fetal human brain tissue, mouse brain tissue, rat brain tissue

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28175-1-AP (APH1B antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.