

For Research Use Only

APH1B Polyclonal antibody, PBS Only

Catalog Number: 28175-1-PBS



Basic Information

Catalog Number:

28175-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14066

GenBank Accession Number:

BC020905

GeneID (NCBI):

83464

UNIPROT ID:

Q8WW43

Full Name:

anterior pharynx defective 1 homolog B (C. elegans)

Calculated MW:

257 aa, 28 kDa

Observed MW:

30 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28175-1-AP (APH1B antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 28175-1-PBS in a different storage buffer formulation.