

For Research Use Only

KTI12 Polyclonal antibody

Catalog Number: 28237-1-AP



Basic Information

Catalog Number: 28237-1-AP	GenBank Accession Number: BC012173	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 400 µg/ml by Nanodrop and 300 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 112970	Recommended Dilutions: WB 1:2000-1:8000
Source: Rabbit	Full Name: KTI12 homolog, chromatin associated (S. cerevisiae)	
Isotype: IgG	Calculated MW: 354 aa, 39 kDa	
Immunogen Catalog Number: AG28201	Observed MW: 39 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, HeLa cells, HepG2 cells, Jurkat cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28237-1-AP (KT112 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.