

For Research Use Only

RAB43 Polyclonal antibody

Catalog Number: 28240-1-AP



Basic Information

Catalog Number: 28240-1-AP	GenBank Accession Number: BC062319	Purification Method: Antigen affinity purification
Size: 150ul , Concentration: 600 µg/ml by Nanodrop and 333 µg/ml by Bradford method using BSA as the standard;	GeneID (NCBI): 339122	Recommended Dilutions: WB 1:500-1:2000
Source: Rabbit	Full Name: RAB43, member RAS oncogene family	
Isotype: IgG	Calculated MW: 212 aa, 23 kDa	
Immunogen Catalog Number: AG26911	Observed MW: 23 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HEK-293 cells, HepG2 cells
Species Specificity: Human	

Background Information

Storage


Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28240-1-AP (RAB43 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.