

For Research Use Only

# Phospho-RIPK3 (Thr231/Ser232) Polyclonal antibody



Catalog Number: 28253-1-AP

## Basic Information

**Catalog Number:**

28253-1-AP

**Size:**

100ul , Concentration: 200 µg/ml by Bradford method using BSA as the standard;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_019955

**GeneID (NCBI):**

56532

**UNIPROT ID:**

Q9QZL0

**Full Name:**

receptor-interacting serine-threonine kinase 3

**Calculated MW:**

53 kDa

**Purification Method:**

Antigen affinity purification

## Applications

**Tested Applications:**

ELISA

**Species Specificity:**

Human

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data